

penreco®

Red LED Booster Serum

Use this innovative serum to boost the effects of red LED light equipment and also to soothe sensitive skin



Red LED Booster Serum

PRODUCTS

Raw Material	Composition	% w/w	Supplier
PHASE A			
Lubrajel™ MS Free	Aqua, Glycerin, Glyceryl Acrylate/Acrylic Acid Copolymer, Phenoxyethanol	41.00	Guardian
Lexgard® O-EZ	Caprylyl Glycol, Propanediol	1.00	Inolex
Versagel® MD500T	Isododecane, Ethylene/Propylene/Styrene Copolymer, Butylene/Ethylene/Styrene Copolymer, Pentaerythryl Tetra-di-t-butyl Hydroxyhydrocinnamate	41.50	Penreco
MultiEx BSASM™	Butylene Glycol, Aqua, Centella Asiatica Extract, Polygonum Cuspidatum Root Extract, Scutellaria Baicalensis Root Extract, Camellia Sinensis Leaf Extract, Glycyrrhiza Glabra Root Extract, Chamomilla Recutita Flower Extract, Rosmarinus Officinalis Leaf Extract	2.50	Clariant AI
Sebocure™	Aqua, Phytosterols, Caprylic/Capric Triglyceride, Lecithin, Squalane, Glycyrrhiza Glabra Root Extract, Soy Isoflavones, Butyrospermum Parkii Butter, Dioscorea Villosa Root Extract, Niacinamide, Pyridoxine HCl, Azelaic Acid, Ceramide 3, Zinc PCA	2.00	Clariant AI
Homeostatine	Aqua, Glycerin, Enteromorpha Compressa Extract, Caesalpina Spinosa Gum, Phenoxyethanol, Potassium Sorbate	2.00	Provital
Isopentylidol	Isopentylidol	10.00	Kuraray

MIXING INSTRUCTIONS

Mix phase A, in order until homogenous between each addition.